





KMedic® Super Cut Scissors



KM 35-450 KM 35-452 KM 35-454 KM 35-460

straight
 KM 35-450 3½" delicate
 KM 35-454 4½"
curved
 KM 35-452 3½" delicate
 KM 35-460 4½"



SUPER-CUT IRIS SCISSORS sharp points



Delicate short narrow tips for dissecting fine tissue.



SUPER-CUT STEVENS SCISSORS

blunt points
straight
 KM 35-484 4½"
curved
 KM 35-490 4½"
 KM 35-494 5¼"





KM 35-220 KM 35-250

SUPER-CUT OPERATING SCISSORS
 KM 35-220 KM 35-250
 5½" straight 5½" curved
 sharp/blunt points sharp/blunt points





Very delicate blunt scissors with long handles for fine tissue dissection.

SUPER-CUT JAMISON SCISSORS
 KM 35-496
 6¼" curved blunt points, very delicate



KM 35-364 KM 35-370


SUPER-CUT BABY METZENBAUM SCISSORS
 KM 35-364 KM 35-370
 4½" straight 4½" curved
 blunt points blunt points




KM 35-382




KM 35-394




KM 35-396



KM 35-406



KM 35-418





KM 35-430


straight
KM 35-380 6"
KM 35-392 7"

curved
KM 35-382 6"
KM 35-394 7"
KM 35-396 7" delicate
KM 35-406 8"
KM 35-418 9"
KM 35-430 11"


SUPER-CUT METZENBAUM SCISSORS blunt points




KM 35-280



KM 35-292



KM 35-316



KM 35-322


SUPER-CUT MAYO SCISSORS

beveled blades, blunt points


<i>straight</i>	<i>curved</i>
KM 35-280 5½"	KM 35-286 5½"
KM 35-292 6¾"	KM 35-298 6¾"
KM 35-316 9"	KM 35-322 9"
	KM 35-340 11"

Super-Cut Bandage Scissors

KMedic developed these remarkably sharp scissors in response to the demands of doctors, nurses, and other O.R. personnel.



KM 35-210




KM 35-212

SUPER-CUT LISTER BANDAGE SCISSORS

KM 35-210
5 1/2"


KM 35-212
7 1/4"



The delicate offset shanks and the sharp Super-Cut edge provide better leverage. The fine knob slides easily under finger and toe bandages.

SUPER-CUT HI-LEVEL BANDAGE SCISSORS

KM 35-202
5 1/2" shanks angled on side



All the benefits of our Super-Cut Lister scissors plus an oversized finger ring, which helps relieve pressure on the third metacarpal.

SUPER-CUT LISTER BANDAGE SCISSORS WITH LARGE FINGER RING

KM 35-214
8"



KM 35-646

TC METZENBAUM SCISSORS blunt points

straight

KM 35-640 6"

KM 35-644 7"

KM 35-652 9"

curved

KM 35-642 6"

KM 35-646 7"

KM 35-654 9"



KM 35-620



KM 35-622

TC MAYO SCISSORS beveled blades

straight

KM 35-620 5 1/2"

KM 35-624 6 3/4"

KM 35-628 9"

curved

KM 35-622 5 1/2"

KM 35-626 6 3/4"

KM 35-630 9"

TELEFLEX MEDICAL OEM

5307 95th Avenue, Kenosha, WI 53144 USA

Toll Free: 800 295 8505 (USA only)

Phone: 262 657 2800

E-mail: oeinfo@teleflexmedical.com

Website: www.teleflexmedicoem.com

©2011 Teleflex Incorporated. All rights reserved. Teleflex and KMedic are registered trademarks of Teleflex Incorporated or its affiliates.