


KMedic® Bite and Handle Styles


IDENTIFICATION CHART

SHAFT LENGTH
measured in closed position




BITE STYLES punch jaws, thin footplate


Forward 40°



Up 90°




Down 90°



BITE STYLES punch jaws, regular footplate


Forward 40°



Up 90°



Down 90°



BITE STYLES cup jaws

Straight



Up




Down




PUNCH JAWS

Foot plate



Bite opening



Bite size



HANDLE STYLES



STYLE I
Ring handle
handle height: 4"



STYLE II
Kerrison Handle
handle height: 6 1/4"



STYLE III
Love-Kerrison Handle
(Cherry)
handle height: 6"



STYLE IV
Large Kerrison Handle
handle height: 6 1/4"

All KMedic spinal rongeurs and punches have bite size etched on handles for easy identification.

HANDLE STYLES



STYLE V
Small Kerrison Handle
handle height: 5"



STYLE VI
Ferris-Smith
Kerrison Ring Handle
handle height: 5"



STYLE VII
Spurling Kerrison Handle
handle height: 4 3/4"
Wider grip improves
comfort and control.



STYLE VIII
Medium Kerrison Handle
handle height: 5 1/4"

All KMedic spinal rongeurs and punches have bite size etched on handles for easy identification.