




















KMedic® Thin Footplate Kerrison Punch



- Thin footplate designed specifically for small bone in the cervical area.
- Gold hexagonal screw denotes thin footplate.
- Bite size etched on handle for quick identification.
- Care should be taken when using delicate thin footplate punches.
- Do not twist or pry instrument while cutting.
- Do not use on large bones.

(Handle style VII)

 Forward 40° Thin Footplate	 Up 90° Thin Footplate	 Down 90° Thin Footplate
 1mm bite forward KM 55-891 6" KM 55-904 7" KM 55-925 8"	 1mm bite up KM 55-905 7"	 1mm bite down KM 55-906 7"
 2mm bite forward KM 55-892 6" KM 55-908 7" KM 55-926 8"	 2mm bite up KM 55-909 7"	 2mm bite down KM 55-910 7"
 3mm bite forward KM 55-895 6" KM 55-911 7" KM 55-927 8"	 3mm bite up KM 55-912 7"	 3mm bite down KM 55-913 7"
 4mm bite forward KM 55-898 6" KM 55-914 7" KM 55-928 8"	 4mm bite up KM 55-915 7"	 4mm bite down KM 55-916 7"
 5mm bite forward KM 55-901 6" KM 55-917 7" KM 55-929 8"	 5mm bite up KM 55-918 7"	 5mm bite down KM 55-919 7"