

BLACK Force Fiber® Suture

Introducing
a Suture Braided
from BLACK
Ultra-High
Molecular Weight
Polyethylene
(UHMWPE) Fiber

- Distinct Color
- High Strength
- Durable
- Low Knot Profile

Sizes (USP)	0	1	2	Additional sizes evaluated upon request
Solid				
• Black	■	■	■	Black-dyed UHMWPE
Cobraids				Additional cobraids evaluated upon request
• White/Black	■	■	■	Cobraid of undyed and dyed UHMWPE
• Blue/Black	■	■	■	Cobraid of dyed UHMWPE
• Black/White	■	■	■	Cobraid of undyed and dyed UHMWPE
Patterns				
• Criss-Cross	■	■	■	
• Striped	■	■	■	
Configurations				
• Spools, 500 Yards	■	■	■	Non-sterile
• Cut Lengths, 12" to 60" (31 to 152 cm)	■	■	■	Non-sterile

Please note: **Black** Force Fiber® Suture is currently only available as a *non-sterile* product.

Teleflex®
MEDICAL OEM

Work With The Suture Experts™

Phone: 1.800.474.0178 or 1.508.964.6021

Email: oeminfo@teleflex.com

Web: www.teleflexmedicaloem.com

©2018, Teleflex Incorporated. All rights reserved. Teleflex, Force Fiber, and "Work With The Suture Experts" are trademarks or registered trademarks of Teleflex Incorporated in the U.S. and/or other countries. 4390-1018