

# Force Fiber® OrthoTape™ Braid

Sizes	1mm	1.5mm	2mm	2.5mm	3.5mm	4mm	5mm	<i>This ultra high-strength tape features a low profile, broad footprint, and high lubricity. It can address the concern of tissue pull-through that sometimes occurs with round suture.</i>
Average Width	0.70 - 1.10mm	1.30 - 1.60mm	1.70 - 2.00mm	2.10 - 2.70mm	3.20 - 3.60mm	3.40 - 4.10mm	4.50 - 5.40mm	
Average Height	0.08 - 0.45mm	0.12 - 0.52mm	0.26 - 0.60mm	0.10 - 1.00mm	0.23 - 1.12mm	0.40 - 1.40mm	0.40 - 1.80mm	
Average Min. Knot Pull	≥ 10 lbs	≥ 14 lbs	≥ 20 lbs	≥ 28 lbs	≥ 34 lbs	≥ 120 lbs*	≥ 150 lbs*	
<b>Solids</b>								
White	■	■	■	■	■	■	■	Undyed UHMWPE
Blue	■	■	■					Blue-dyed UHMWPE
<b>Multi-Colored</b>								
White/Blue		■						Cobraid of undyed and blue UHMWPE
White/Green		■	■	■	■			Cobraid of undyed UHMWPE and green polyethylene terephthal multifilament (PET)
<b>Patterns</b>								
Criss-Cross or Tracer Strands		■	■	■	■			Custom patterns evaluated upon request
<b>Configurations</b>								
Spools, 250 yards	■	■	■	■	■	■	■	Non-sterile
Cut Length, Tipped	■	■	■	■				<ul style="list-style-type: none"> <li>• Non-sterile or sterile</li> <li>• With or without needles</li> </ul>
Cut Length, Heat Stiffened	■	■	■					<ul style="list-style-type: none"> <li>• Non-sterile or sterile</li> <li>• With or without needles</li> </ul>

\*Please note: Pass/fail test