



Safety OTN
Peelable Introducers

TFX OEM® OVER-THE-NEEDLE PEELABLE INTRODUCERS

CUSTOM OTN INTRODUCERS ARE AVAILABLE IN TWO STYLES

- Safety over-the-needle (OTN) introducer has a passive safety device that is activated automatically.
- Both safety and standard OTN introducers are available in sizes 2F to 6F, (approximately 22 GA to 14 GA needle size). *See chart on opposite side.*
- Applications: Percutaneous insertion of small PICC, midline, or other peripheral catheters and pacing leads.
- Clear, large volume, flashback chamber allows for visualization and blood containment.
- Flashback plug permits use of syringe.
- Funneled hub design for smooth insertion of PICC, midline, or other catheters.
- Peel-apart introducer for easy removal.
- Ribbed wing design for secure handling.
- Needle is silicone coated to minimize resistance during insertion.
- Smooth transition between needle and the sheath to minimize tissue trauma.
- Cover tube protects the needle prior to use.



Standard OTN Peelable Introducers



PARTNER WITH TELEFLEX MEDICAL OEM

Whether you need a component or a finished medical device, Teleflex Medical OEM should be your product design, development, and production partner. Along with our “can-do” attitude, we offer engineering expertise, extensive prototyping and manufacturing capabilities, regulatory services, and a dedication to customer service.

SAFETY OTN PEELABLE INTRODUCERS

- Personnel are protected by a passive safety device that is activated automatically.
- Meets OSHA safety requirements.
- Sizes 2F to 6F, (approximately 22 GA to 14 GA needle size).
- Custom configuration and labeling are available.

STANDARD OTN PEELABLE INTRODUCERS

- Sizes 2F to 6F, (approximately 22 GA to 14 GA needle size).
- Custom configuration and labeling are available.

Dimensions

Needle Size	Introducer Size	I.D./Sheath Introducer (Min)	O.D./Sheath Introducer
22 GA	2F	0.030" (0.76mm)	0.044" (1.12mm)
19 GA	3F	0.044" (1.12mm)	0.058" (1.47mm)
17 GA	4F	0.060" (1.52mm)	0.074" (1.88mm)
15 GA	5F	0.074" (1.88mm)	0.088" (2.24mm)
14 GA	6F	0.086" (2.18mm)	0.100" (2.54mm)

TELEFLEX MEDICAL OEM

50 Plantation Drive
Jaffrey, NH 03452 USA

Toll free: +1.800.548.6600
Phone: +1.603.532.7706

E-mail:
oeminfo@teleflex.com

Website:
www.teleflexmedialoem.com

Imagine it. Achieve it.®

©2012 Teleflex Incorporated.
All rights reserved.
Teleflex, TFX OEM, and “Imagine It. Achieve it.” are registered trademarks of Teleflex Incorporated.
1047-0612