

BLACK Force Fiber® Suture

- Braided from 100% Ultra-High Molecular Weight Polyethylene (UHMWPE)
- Distinct Color
- Durable
- Low Knot Profile
- High Strength
- This Suture can be Custom Configured to Meet the Requirements of Your Application, or Incorporated into an Anchor or Driver of Your Design.

Indications

Sterile Force Fiber® Suture is 510k cleared as a Class II device in the USA for approximation and/or ligation of soft tissue, including the use in cardiovascular surgeries and use of allograft tissue for orthopedic surgeries.

Sterile Force Fiber® Suture is also CE marked in the EU as a class IIb device, and licensed in Canada as a class III device, for approximation and/or ligation of soft tissue including the use of allograft tissue for orthopedic surgeries.

Sterile Force Fiber® Suture is contraindicated in the European Union and Canada for cardiovascular use.

Sizes (USP)	0	1	2	Additional sizes evaluated upon request
Solid				
• Black	■	■	■	Dyed 100% UHMWPE
Multi-Colored				
<i>Additional cobraids evaluated upon request</i>				
• White/Black	■	■	■	Undyed and dyed 100% UHMWPE
• Blue/Black	■	■	■	Dyed 100% UHMWPE
• Black/White	■	■	■	Dyed and undyed 100% UHMWPE
• White/Blue/Black	■	■	■	Undyed and dyed 100% UHMWPE
Patterns				
• Crisscross	■	■	■	
• Striped	■	■	■	
Configurations				
• Spools, 500 yards	■	■	■	Non-sterile only
• Cut lengths, 12" to 60" (31 to 152 cm)	■	■	■	Sterile, private label; non-sterile <i>Additional cut lengths evaluated upon request</i>
• Cut lengths, with or without tipping	■	■	■	Sterile, private label; bulk non-sterile
• Cut lengths, with attached needles	■	■	■	Sterile, private label; bulk non-sterile
• Preloaded in assembled medical devices	■	■	■	

Teleflex®
MEDICAL OEM

Phone: 1.800.474.0178 or 1.508.964.6021

Email: oeinfo@teleflex.com

Web: www.teleflexmedicaloem.com