



# QUICK-TURN OTW SEMI-CUSTOM BALLOON CATHETERS



ITEM NO.	DESCRIPTION
1	5Fr Co-Axial Shaft Pebax® 72D, 295C Cobalt Blue
2	Y-Arm Hub, Polycarbonate ISO 80369 Compliant*
3	Strain Relief, TPE
4	Balloon, Nylon 12/Vestamid ML 21
5	RO Marker Band(s), Pt 90% / Ir 10%

Balloon Size	†
Diameter (mm)	3, 4, 5, 6, 7
Body Length (mm)	20, 40, 80, 100
Catheter	
Working Length (cm)	0-250
Distal Tip Length (mm)	5-15
Strain Relief Color	295 C Cobalt Blue or 6C Cool Grey
Marker Band(s) options	0 1 (Centered) 2 (Each Cone Angle)
Additional Options	
Balloon Wrapped	
Balloon Protector	






**LET'S GET TO WORK  
ON YOUR PROJECT**

- Phone: 1.800.474.0178 (USA only) or 1.508.964.6021
- E-mail: [oeminfo@teleflex.com](mailto:oeminfo@teleflex.com)
- Web: [www.teleflexmedicaloem.com](http://www.teleflexmedicaloem.com)

*Work With The Experts™*

**Teleflex®**  
MEDICAL OEM

[teleflexmedicaloem.com](http://teleflexmedicaloem.com)   

\*ISO 80369 Luer testing on file with Teleflex Medical OEM.  
† For additional custom balloon sizes and lengths not listed, please contact your local sales representative for quote.

Customer acknowledges that Teleflex is providing the above product in support of defining customers product design requirements.

Responsibility of the Customer:

- Ownership of product design, biocompatibility testing, and risk management.
- Regulatory responsibility and ownership associated with this product.
- Export controls and requirements, if applicable.

Teleflex, the Teleflex logo, and Work With The Experts are trademarks or registered trademarks of Teleflex Incorporated or its affiliates in the U.S. and/or other countries. Other names are trademarks of their respective owners.